3M™ Clean-Trace™ ATP Monitoring System

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- Surfaces are disinfected or sterilized in order to prevent transmission of pathogens.
- All soil must be removed from surfaces as most disinfectants and sterilants do not penetrate through organic matter.
- Soil can protect microbes from the action of disinfection and sterilization.
- Surviving microbes have potential to be transmitted to another patient causing infection.
- Recommended Practices, Standards and Guidelines are beginning to emphasize the importance of monitoring cleaning efficacy.
 - FDA, CDC, AORN, APIC, SHEA, AAMI

Know the difference

Cleaning

- Removal of organic soil
- Microbes and soil can still be present
- Device can still be infectious

High-Level Disinfection (HLD)

- Microbial kill under defined conditions
- Spores are not killed
- Effectiveness dependent on meticulous cleaning

Sterilization

- Kills all living organisms including spores
- Effectiveness dependent on meticulous cleaning

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
The background is a complex, abstract pattern of overlapping triangles and polygons in various shades of green, ranging from light lime green to dark forest green. The shapes are irregular and create a sense of depth and movement.

3M™ Clean-Trace™ ATP Monitoring System

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3M™ Clean-Trace™ Systém komponenty k použití

Hardware
3M™ Clean-Trace
Luminometer LX25

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3M™ ATP Surface Test UXC	3M™ ATP Water Test H2O
-----------------------------	---------------------------

Software

3M™ Quality Control Data Manager (QCDM)

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
Quality Control Data
Manager Dashboard
Examples

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What is one of the biggest risks to a patient after being admitted to a hospital?

There is a 50% chance that you will catch a multi-drug resistant organism (MDRO) from your patient room if it was previously occupied by a patient who was colonized or infected with a multi-drug resistant organism (MDRO).

— Otter et al. *Am J Infect Control* 2013;41(5 Suppl):S6-11.

Factors that contribute to patient risk


- Transmission of MDRO from the environment to patients
- Inadequate cleaning
- The environment is now harder to clean than in the past
- Multidrug Resistant Organisms (MDRO) and their persistence in the environment

Clostridium difficile
>5 months to years

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

Monitoring of High Risk Areas and Equipment

Pass/Fail Threshold

- ≤ 250 RLU = Pass
- ≥ 251 = Fail


Test point plan is available

Frequency

- Monitor after every cleaning, every time.

Chirurgické nástroje

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Pass/Fail Threshold

*After Instrument Removed from
Washer-Disinfector*

Pass ≤ 150 RLU

Fail ≥ 151 RLU

Rutinni Monitoring

- ✓ Monitor 10 hard to clean instruments
- ✓ Monitor every day, every shift
- ✓ Monitor after Washer-Disinfector Step
- ✓ Pass ≤ 150 RLU Fail ≥ 151 RLU
- ✓ If more than 2 instruments fail, and they are in the same set, the entire set should be re-cleaned
- ✓ If instruments are from different sets, then tracking and trending over time is essential to catch any issues

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- Complex design
- Multiple, long, narrow channels that are difficult to clean
- Lack of consistent effective training
- Lack of time and resources for adequate reprocessing
- Visual inspection not adequate to monitor efficacy of reprocessing
- >130 steps involved in reprocessing!

Kde v procesu 3M™ Clean-Trace™ používá?

After manual Cleaning

Before High-Level Disinfection

Not recommended after high High-Level Disinfection or Sterilization or after Drying/Storage

Pass/Fail Threshold

Pass ≤ 200 RLU

Fail ≥ 201 RLU

- Used for all test points, all scope types.
- Only threshold for ATP monitoring that is referenced in AAMI ST91*

Recommended test points—can vary by endoscope type

For ALL Scope Types

1. The outer distal end (ATP Surface Test UXC)
2. The suction biopsy channel (ATP Water Test H2O)

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AAMI (Association for the Advancement for Medical Instrumentation)

ANSI/AAMI ST91:2015. Flexible and semi-rigid endoscope processing in health-care facilities.

AORN (Association of periOperative Registered Nurses)

Guidelines for Processing Flexible Endoscopes

SGNA (Society for Gastroenterology Nurses and Associates)

Standards of Infection Prevention in Reprocessing Flexible Gastrointestinal Endoscopes

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3M™ Clean-Trace™ ATP Monitoring System zahrnuje

- 3M™ Clean-Trace™ Luminometer LX25
- 3M™ Quality Control Data Manager
- 3M™ Clean-Trace™ Surface Test UXC
- 3M™ Clean-Trace™ Water Test H2O

Thank you

lsyrovatkova@mmm.com